



scn Antibody

Product Code	CSB-PA858805HA01SKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99SU9
Immunogen	Recombinant Staphylococcus aureus Staphylococcal complement inhibitor protein (32-116AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Staphylococcal complement inhibitor (SCIN), scn
Immunogen Species	Staphylococcus aureus (strain N315)
Research Area	Others
Target Names	scn
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.