



scn Antibody, Biotin conjugated

Product Code	CSB-PA858805HD01SKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99SU9
Immunogen	Recombinant Staphylococcus aureus Staphylococcal complement inhibitor protein (32-116AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Staphylococcal complement inhibitor (SCIN), scn
Immunogen Species	Staphylococcus aureus (strain N315)
Research Area	Others
Target Names	scn