



PAIP2 Antibody, Biotin conjugated

Product Code	CSB-PA859027LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPZ3
Immunogen	Recombinant Human Polyadenylate-binding protein-interacting protein 2 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyadenylate-binding protein-interacting protein 2 (PABP-interacting protein 2) (PAIP-2) (Poly(A)-binding protein-interacting protein 2), PAIP2, PAIP2A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PAIP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.