



PAIP2 Antibody, Biotin conjugated

Product Code	CSB-PA859027LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPZ3
Immunogen	Recombinant Human Polyadenylate-binding protein-interacting protein 2 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyadenylate-binding protein-interacting protein 2 (PABP-interacting protein 2) (PAIP-2) (Poly(A)-binding protein-interacting protein 2), PAIP2, PAIP2A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PAIP2