



SPIN2A Antibody

Product Code	CSB-PA859525LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99865
Immunogen	Recombinant Human Spindlin-2A protein (1-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spindlin-2A (Protein DXF34) (Spindlin-like protein 2A) (SPIN-2) (SPIN-2A), SPIN2A, DXF34 SPIN2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPIN2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.