



SPIN2A Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA859525LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99865 |
| Immunogen | Recombinant Human Spindlin-2A protein (1-258AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Spindlin-2A (Protein DXF34) (Spindlin-like protein 2A) (SPIN-2) (SPIN-2A), SPIN2A, DXF34 SPIN2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | SPIN2A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |