



S1PR3 Antibody, HRP conjugated

Product Code	CSB-PA859932LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99500
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 3 protein (303-378AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingosine 1-phosphate receptor 3 (S1P receptor 3) (S1P3) (Endothelial differentiation G-protein coupled receptor 3) (Sphingosine 1-phosphate receptor Edg-3) (S1P receptor Edg-3), S1PR3, EDG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	S1PR3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.