



# S1PR3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA859932LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q99500
<b>Immunogen</b>	Recombinant Human Sphingosine 1-phosphate receptor 3 protein (303-378AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Sphingosine 1-phosphate receptor 3 (S1P receptor 3) (S1P3) (Endothelial differentiation G-protein coupled receptor 3) (Sphingosine 1-phosphate receptor Edg-3) (S1P receptor Edg-3), S1PR3, EDG3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	S1PR3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.