



TSC22D3 Antibody, Biotin conjugated

Product Code	CSB-PA860343LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99576
Immunogen	Recombinant Human TSC22 domain family protein 3 protein (58-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TSC22 domain family protein 3 (DSIP-immunoreactive peptide) (Protein DIP) (hDIP) (Delta sleep-inducing peptide immunoreactor) (Glucocorticoid-induced leucine zipper protein) (GILZ) (TSC-22-like protein) (TSC-22-related protein) (TSC-22R), TSC22D3, DSIPI GILZ
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TSC22D3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.