



P2RY13 Antibody, HRP conjugated

Product Code	CSB-PA860651LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPV8
Immunogen	Recombinant Human P2Y purinoceptor 13 protein (162-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2Y purinoceptor 13 (P2Y13) (G-protein coupled receptor 86) (G-protein coupled receptor 94), P2RY13, GPR86 GPR94
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RY13