



P2RY13 Antibody, Biotin conjugated

Product Code	CSB-PA860651LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPV8
Immunogen	Recombinant Human P2Y purinoceptor 13 protein (162-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2Y purinoceptor 13 (P2Y13) (G-protein coupled receptor 86) (G-protein coupled receptor 94), P2RY13, GPR86 GPR94
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RY13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.