



RHOXF2 Antibody, HRP conjugated

Product Code	CSB-PA861117LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQY4
Immunogen	Recombinant Human Rhox homeobox family member 2 protein (1-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhox homeobox family member 2 (Paired-like homeobox protein PEPP-2) (Testis homeobox gene 1), RHOXF2, PEPP2 THG1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RHOXF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.