



PSMG3 Antibody, HRP conjugated

Product Code	CSB-PA861132LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BT73
Immunogen	Recombinant Human Proteasome assembly chaperone 3 protein (1-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome assembly chaperone 3 (PAC-3) (hPAC3), PSMG3, C7orf48 PAC3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMG3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.