



ATG10 Antibody, Biotin conjugated

Product Code	CSB-PA861988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0Y0
Immunogen	Recombinant Human Ubiquitin-like-conjugating enzyme ATG10 protein (1-50AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like-conjugating enzyme ATG10 (EC 2.3.2.-) (Autophagy-related protein 10) (APG10-like), ATG10, APG10L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ATG10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.