



KCNK13 Antibody, HRP conjugated

Product Code	CSB-PA862058LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB14
Immunogen	Recombinant Human Potassium channel subfamily K member 13 protein (292-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 13 (Tandem pore domain halothane-inhibited potassium channel 1) (THIK-1), KCNK13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNK13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.