



SCAPER Antibody, Biotin conjugated

Product Code	CSB-PA863150LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY12
Immunogen	Recombinant Human S phase cyclin A-associated protein in the endoplasmic reticulum protein (1109-1400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S phase cyclin A-associated protein in the endoplasmic reticulum (S phase cyclin A-associated protein in the ER) (Zinc finger protein 291), SCAPER, KIAA1454 ZNF291
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCAPER
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.