



KRTAP3-1 Antibody, Biotin conjugated

Product Code	CSB-PA863160LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYR8
Immunogen	Recombinant Human Keratin-associated protein 3-1 protein (2-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 3-1 (High sulfur keratin-associated protein 3.1) (Keratin-associated protein 3.1), KRTAP3-1, KAP3.1 KRTAP3.1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRTAP3-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.