



PPIL3 Antibody, FITC conjugated

Product Code	CSB-PA863962LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2H8
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase-like 3 protein (2-161AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase-like 3 (PPIase) (EC 5.2.1.8) (Cyclophilin J) (CyPJ) (Cyclophilin-like protein PPIL3) (Rotamase PPIL3), PPIL3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPIL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.