





ADAM7 Antibody

Product Code	CSB-PA863964LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2U9
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 7 protein (177-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 7 (ADAM 7) (Sperm maturation-related glycoprotein GP-83), ADAM7, GP83
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ADAM7

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA863964LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA863964LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.