



MFSD1 Antibody, HRP conjugated

Product Code	CSB-PA863974LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3U5
Immunogen	Recombinant Human Major facilitator superfamily domain-containing protein 1 protein (1-38AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major facilitator superfamily domain-containing protein 1 (Smooth muscle cell-associated protein 4) (SMAP-4), MFSD1, SMAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MFSD1