



POMK Antibody, HRP conjugated

Product Code	CSB-PA863985LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5K3
Immunogen	Recombinant Human Protein O-mannose kinase protein (44-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein O-mannose kinase (POMK) (EC 2.7.1.183) (Protein kinase-like protein SgK196) (Sugen kinase 196), POMK, SGK196
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POMK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.