



POMK Antibody, HRP conjugated

Product Code	CSB-PA863985LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5K3
Immunogen	Recombinant Human Protein O-mannose kinase protein (44-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein O-mannose kinase (POMK) (EC 2.7.1.183) (Protein kinase-like protein SgK196) (Sugen kinase 196), POMK, SGK196
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POMK