





ASB7 Antibody

Product Code	CSB-PA863988LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H672
Immunogen	Recombinant Human Ankyrin repeat and SOCS box protein 7 protein (1-274AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin repeat and SOCS box protein 7 (ASB-7), ASB7
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ASB7

Image



Western blot
All lanes: ASB7 antibody at 8µg/ml + A431 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 37, 31 kDa
Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA863988LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.