





ASB7 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA863988LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H672 |
| Immunogen | Recombinant Human Ankyrin repeat and SOCS box protein 7 protein (1-274AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ankyrin repeat and SOCS box protein 7 (ASB-7), ASB7 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ASB7 |

Image



Western blot
All lanes: ASB7 antibody at 8µg/ml + A431 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 37, 31 kDa
Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA863988LA01HU at dilution of 1:100

