



CIAO3 Antibody, Biotin conjugated

Product Code	CSB-PA863991LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6Q4
Immunogen	Recombinant Human Cytosolic iron-sulfur assembly component 3 protein (2-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic iron-sulfur assembly component 3, Cytosolic Fe-S cluster assembly factor NARFL (Iron-only hydrogenase-like protein 1) (IOP1) (Nuclear prelamina A recognition factor-like protein) (Protein related to Narf), CIAO3, NARFL, PRN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CIAO3