



GAREM1 Antibody, Biotin conjugated

Product Code	CSB-PA863993LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H706
Immunogen	Recombinant Human GRB2-associated and regulator of MAPK protein (12-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRB2-associated and regulator of MAPK protein 1 (GRB2-associated and regulator of MAPK1), GAREM1, C18orf11 FAM59A GAREM
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GAREM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.