



BCO1 Antibody, Biotin conjugated

Product Code	CSB-PA864015LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAY6
Immunogen	Recombinant Human Beta,beta-carotene 15,15'-dioxygenase protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta,beta-carotene 15,15'-dioxygenase (EC 1.13.11.63) (Beta-carotene dioxygenase 1) (Beta-carotene oxygenase 1), BCO1, BCDO BCDO1 BCMO1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BCO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.