




CABP5 Antibody

Product Code	CSB-PA865092LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP86
Immunogen	Recombinant Human Calcium-binding protein 5 protein (1-173AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding protein 5 (CaBP5), CABP5, CABP3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CABP5

Image



Western Blot

All lanes: CABP5 antibody at 12µg/ml

Lane 1: Mouse kidney tissue

Lane 2: Mouse thymus tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 20 kDa

Observed band size: 20 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.