



GPR84 Antibody, FITC conjugated

Product Code	CSB-PA865107LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQS5
Immunogen	Recombinant Human G-protein coupled receptor 84 protein (202-320AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor 84 (Inflammation-related G-protein coupled receptor EX33), GPR84, EX33
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR84
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.