



GJB4 Antibody, HRP conjugated

Product Code	CSB-PA865136LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTQ9
Immunogen	Recombinant Human Gap junction beta-4 protein (211-266AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gap junction beta-4 protein (Connexin-30.3) (Cx30.3), GJB4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GJB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.