



SPTLC3 Antibody, Biotin conjugated

Product Code	CSB-PA865142LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUV7
Immunogen	Recombinant Human Serine palmitoyltransferase 3 protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine palmitoyltransferase 3 (EC 2.3.1.50) (Long chain base biosynthesis protein 2b) (LCB2b) (Long chain base biosynthesis protein 3) (LCB 3) (Serine-palmitoyl-CoA transferase 3) (SPT 3), SPTLC3, C20orf38 SPTLC2L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPTLC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.