



SPTLC3 Antibody, Biotin conjugated

Product Code	CSB-PA865142LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUV7
Immunogen	Recombinant Human Serine palmitoyltransferase 3 protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine palmitoyltransferase 3 (EC 2.3.1.50) (Long chain base biosynthesis protein 2b) (LCB2b) (Long chain base biosynthesis protein 3) (LCB 3) (Serine-palmitoyl-CoA transferase 3) (SPT 3), SPTLC3, C20orf38 SPTLC2L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPTLC3