




FKBP14 Antibody

Product Code	CSB-PA865159LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWM8
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP14 protein (20-211AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase FKBP14 (PPIase FKBP14) (EC 5.2.1.8) (22 kDa FK506-binding protein) (22 kDa FKBP) (FKBP-22) (FK506-binding protein 14) (FKBP-14) (Rotamase), FKBP14, FKBP22
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FKBP14

Image



Western blot

All lanes: FKBP14 antibody at 1µg/ml

Lane 1: Mouse stomach tissue

Lane 2: Mouse gonadal tissue

Lane 3: U937 whole cell lysate

Lane 4: PC-3 whole cell lysate

Lane 5: 293T whole cell lysate

Lane 6: THP-1 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.