



FKBP14 Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA865159LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NWM8 |
| Immunogen | Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP14 protein (20-211AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidyl-prolyl cis-trans isomerase FKBP14 (PPIase FKBP14) (EC 5.2.1.8) (22 kDa FK506-binding protein) (22 kDa FKBP) (FKBP-22) (FK506-binding protein 14) (FKBP-14) (Rotamase), FKBP14, FKBP22 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | FKBP14 |