



FKBP14 Antibody, Biotin conjugated

Product Code	CSB-PA865159LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWM8
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP14 protein (20-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase FKBP14 (PPIase FKBP14) (EC 5.2.1.8) (22 kDa FK506-binding protein) (22 kDa FKBP) (FKBP-22) (FK506-binding protein 14) (FKBP-14) (Rotamase), FKBP14, FKBP22
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FKBP14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.