



HCFC1R1 Antibody, Biotin conjugated

Product Code	CSB-PA865161LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWW0
Immunogen	Recombinant Human Host cell factor C1 regulator 1 protein (1-119AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Host cell factor C1 regulator 1 (HCF-1 beta-propeller-interacting protein), HCFC1R1, HPIP
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	HCFC1R1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.