



OSER1 Antibody, Biotin conjugated

Product Code	CSB-PA865164LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX31
Immunogen	Recombinant Human Oxidative stress-responsive serine-rich protein 1 protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxidative stress-responsive serine-rich protein 1 (Oxidative stress-responsive protein 1) (Peroxide-inducible transcript 1 protein), OSER1, C20orf111
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OSER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.