



GPR83 Antibody, FITC conjugated

Product Code	CSB-PA865180LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYM4
Immunogen	Recombinant Human Probable G-protein coupled receptor 83 protein (349-423AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 83 (G-protein coupled receptor 72), GPR83, GPR72 KIAA1540
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPR83
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.