




FBXL17 Antibody

Product Code	CSB-PA866216LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UF56
Immunogen	Recombinant Human F-box/LRR-repeat protein 17 protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/LRR-repeat protein 17 (F-box and leucine-rich repeat protein 17) (F-box only protein 13), FBXL17, FBL17 FBX13 FBXO13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXL17

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA866216LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.