



FBXL17 Antibody, HRP conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA866216LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UF56 |
| Immunogen | Recombinant Human F-box/LRR-repeat protein 17 protein (1-303AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | F-box/LRR-repeat protein 17 (F-box and leucine-rich repeat protein 17) (F-box only protein 13), FBXL17, FBL17 FBX13 FBXO13 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | FBXL17 |