



FBXL17 Antibody, Biotin conjugated

Product Code	CSB-PA866216LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UF56
Immunogen	Recombinant Human F-box/LRR-repeat protein 17 protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/LRR-repeat protein 17 (F-box and leucine-rich repeat protein 17) (F-box only protein 13), FBXL17, FBL17 FBX13 FBXO13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXL17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.