



MIOX Antibody, Biotin conjugated

Product Code	CSB-PA866220LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGB7
Immunogen	Recombinant Human Inositol oxygenase protein (1-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol oxygenase (EC 1.13.99.1) (Aldehyde reductase-like 6) (Kidney-specific protein 32) (Myo-inositol oxygenase) (MI oxygenase) (Renal-specific oxidoreductase), MIOX, ALDRL6 KSP32 RSOR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MIOX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.