



MIOX Antibody, Biotin conjugated

Product Code	CSB-PA866220LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGB7
Immunogen	Recombinant Human Inositol oxygenase protein (1-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inositol oxygenase (EC 1.13.99.1) (Aldehyde reductase-like 6) (Kidney-specific protein 32) (Myo-inositol oxygenase) (MI oxygenase) (Renal-specific oxidoreductase), MIOX, ALDRL6 KSP32 RSOR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MIOX