



NICN1 Antibody, HRP conjugated

Product Code	CSB-PA866273LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSH3
Immunogen	Recombinant Human Nicotin-1 protein (1-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotin-1 (NPCEDRG) (Tubulin polyglutamylase complex subunit 5) (PGs5), NICN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NICN1