




PRRG3 Antibody

Product Code	CSB-PA866327LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZD7
Immunogen	Recombinant Human Transmembrane gamma-carboxyglutamic acid protein 3 protein (102-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane gamma-carboxyglutamic acid protein 3 (Proline-rich gamma-carboxyglutamic acid protein 3) (Proline-rich Gla protein 3), PRRG3, PRGP3 TMG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRRG3

Image



Immunofluorescent analysis of HeLa cells using CSB-PA866327LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)