



PRRG3 Antibody, Biotin conjugated

Product Code	CSB-PA866327LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZD7
Immunogen	Recombinant Human Transmembrane gamma-carboxyglutamic acid protein 3 protein (102-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane gamma-carboxyglutamic acid protein 3 (Proline-rich gamma-carboxyglutamic acid protein 3) (Proline-rich Gla protein 3), PRRG3, PRGP3 TMG3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRRG3