



Gkn2 Antibody

Product Code	CSB-PA866591LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9CQS6
Immunogen	Recombinant Mouse Gastrokine-2 protein (21-184AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrokine-2 (Blottin), Gkn2, Blot
Immunogen Species	Mus musculus (Mouse)
Research Area	Cancer
Target Names	Gkn2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.