



PCED1A Antibody, HRP conjugated

Product Code	CSB-PA867128LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1Q7
Immunogen	Recombinant Human PC-esterase domain-containing protein 1A protein (6-312AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PC-esterase domain-containing protein 1A (Protein FAM113A) (Sarcoma antigen NY-SAR-23), PCED1A, C20orf81 FAM113A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PCED1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.