





TUBB1 Antibody

Product Code	CSB-PA867148LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4B7
Immunogen	Recombinant Human Tubulin beta-1 chain protein (270-451AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin beta-1 chain, TUBB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TUBB1


Image




Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA867148LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA867148LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA867148LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA867148LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.