



TUBB1 Antibody, Biotin conjugated

Product Code	CSB-PA867148LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4B7
Immunogen	Recombinant Human Tubulin beta-1 chain protein (270-451AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin beta-1 chain, TUBB1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TUBB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.